

---

# **xraydb**

***Release 4.4.6***

**Matthew Newville**

**2021-May-03**



# CONTENTS

<b>1</b>	<b>Table of Contents</b>	<b>3</b>
1.1	Installation	3
1.1.1	Development Version	3
1.1.2	Testing	3
1.1.3	Copyright, Licensing, and Re-distribution	4
1.2	X-ray Periodic Table of the Elements	4
1.3	Example Calculations of X-ray properties of materials	5
1.3.1	X-ray attenuation by elements	5
1.3.2	$\mu$ calculations for materials	7
1.3.3	X-ray flux calculations for ionization chambers and photodiodes	10
1.3.4	X-ray mirror reflectivities	12
1.3.5	Darwin widths of monochromator crystals	13
1.4	Using XrayDB from Python	16
1.4.1	The Python xraydb module	16
1.4.2	Atomic Properties	18
1.4.3	Elastic Scattering Factors	18
1.4.4	X-ray Edges	19
1.4.5	X-ray Emission Lines	20
1.4.6	Absorption and Scattering Cross-sections	21
1.4.7	Chemical and Materials database	24
1.4.8	X-ray properties of materials	26
1.5	Overview of Atomic and X-ray Data	32
1.5.1	Elements	32
1.5.2	Physical Units	32
1.5.3	X-ray Edges	33
1.5.4	X-ray Lines	33
1.5.5	Cross Sections	34
1.6	Using the XrayDB xraydb.sqlite	34
1.6.1	Overall Database Schema	34
1.6.2	Version Table	36
1.6.3	Elements Table	36
1.6.4	Xray_Levels Table	36
1.6.5	Xray_Transitions Table	36
1.6.6	Photoabsorption Table	37
1.6.7	Scattering Table	37
1.6.8	Coster_Kronig Table	38
1.6.9	Waasmaier Table	38
1.6.10	KeskiRahkonen_Krause Table	38
1.6.11	Krause_Oliver Table	39
1.6.12	Chantler Table	39

1.7	References . . . . .	39
	<b>Bibliography</b>	<b>41</b>
	<b>Python Module Index</b>	<b>43</b>
	<b>Index</b>	<b>45</b>

XrayDB provides atomic data, characteristic X-ray energies, and X-ray cross sections for the elements in an SQLite3 database, `xraydb.sqlite`. This file can be used directly with SQLite [\[\[Hipp \(2012\)\]\]](#) using standard SQL, or from the many programming language that support SQLite. Some of the components of the database hold arrays of numbers which are stored in the database as JSON-encoded strings, and will need to be unpacked to be used. A Python module providing such an interface is provided [here](#).

The project began with the data from the compilation of basic atomic properties and X-ray absorption edge energies, emission energies, and absorption cross sections from [\[\[Elam, Ravel, and Sieber \(2002\)\]\]](#), who assembled data from a several sources. More data has been added from other sources. Energy widths of core holes for excited electronic levels from [\[\[Keski-Rahkonen and Krause \(1974\)\]\]](#) and [\[\[Krause and Oliver \(1979\)\]\]](#). Elastic X-ray scattering data,  $f_0(q)$  is taken from [\[\[Waasmaier and Kirfel \(1995\)\]\]](#). Resonant scattering cross sections  $f'(E)$  and  $f''(E)$  and absorption cross sections from [\[\[Chantler \(2000\)\]\]](#) as from the [FFAST webpage](#) (but on a finer energy grid, data from [\[\[Chantler \(2016\)\]\]](#)) are also included.

The values here use units that are SI and commonly used in X-ray work. Cross sections are in  $\text{cm}^2/\text{gr}$ , and energies are in eV. Energy-dependent data for cross-sections are typically most reliable between about 250 eV to about 250,000 eV. Elements from  $Z=1$  to 92 are supported, with some data are included for elements between  $Z=93$  and  $Z=98$ .

The current version of the XrayDB is **8.1**, and the version of the Python module is 4.4.6. See Also:

- [XrayDB Web App \(CARS, U Chicago\)](#) and [XrayDB Web App \(xrayabsorption.org\)](#) are interactive web applications to browse the data in this database and make plots of X-ray attenuation, scattering factors, mirror reflectivity, and more.
- Printable Poster-sized *X-ray Periodic Table of the Elements*.
- [XrayDB Github Page](#) has data sources, code, development and issues.
- [PDF Version of this documentation](#)



## TABLE OF CONTENTS

### 1.1 Installation

The X-ray database is held in the SQLite3 file `xraydb.sqlite`. If you are looking for direct use with SQLite, you can download this from [xraydb.sqlite](#).

If you want to use XrayDB from Python, install the XrayDB Python module (which includes the sqlite database), with:

```
pip install xraydb
```

If you are using Anaconda Python, you can also install with:

```
conda install -c gsecars xraydb
```

Depending on your system and Python installation, you may need administrative privileges to install any python library. For many linux and Mac OS X systems, you may need to use *sudo*.

---

**Note:** The Python module supports Python 3.5 and above.

---

#### 1.1.1 Development Version

To work with the data sources or to add or modify data in the XrayDB, you will want to clone or download the full source code kit [xrayDB on github.com](#) which contains the current database, original source data, python module, and files for generating the Periodic Table posters. To get the latest development version, use:

```
git clone https://github.com/xraypy/XrayDB.git
```

#### 1.1.2 Testing

There are a set of tests scripts for the Python interface that can be run with the [pytest](#) testing framework. These are located in the `python/tests` folder. These tests are automatically run as part of the development process. For any release or any master branch from the git repository, running `pytest` should run all of these tests to completion without errors or failures.

## 1.1.3 Copyright, Licensing, and Re-distribution

### Public Domain

To the extent possible, the data, code, and documentation here are placed in the public domain. No claim of copyright is made by the authors and no restrictions of any kind placed on the use of the work here.

The original sources of the data are mostly based on published works with the clear intent of providing data to the general public. Some of the datasets have may not have clear statements of copyright or license, but have been freely available for many years. The work here is a compilation and reformatting of those datasets.

If you use these resources and wish to cite the original data source, please see the documentation for more details.

## 1.2 X-ray Periodic Table of the Elements


XrayDB has been used to generate X-ray Periodic Tables of the Elements.

There is a choice of two sizes, and a choice of an image of one of four prominent scientists associated with the Periodic Table and X-ray spectroscopies.

**Table of X-ray Periodic Tables of the Elements** The poster comes in two sizes: Large is 127.5x61 cm (about 50x24 inches) and Medium is 91x45.5 cm (about 36x17 inches).

Scientist	Large Periodic Table	Medium Periodic Table
Charles Barkla	Barkla (Large)	Barkla (Medium)
Marie Curie	Curie (Large)	Curie (Medium)
Dmitri Mendeleev	Mendeleev (Large)	Mendeleev (Medium)
Henry Moseley	Moseley (Large)	Moseley (Medium)

In general, these table look like this:



but you should definitely download the PDFs linked for high-quality results.



## 1.3 Example Calculations of X-ray properties of materials

Here, a few detailed examples of using the *xraydb.sqlite* to calculate the X-ray properties of materials are shown. These all use the functions in the python *xraydb* module, which is describe in more detail in the next chapter, *Using XrayDB from Python*. The examples will explore some aspects of X-ray physics, but will not give a complete tutorial on the concepts here. For reference see [[Als-Nielsen and McMorrow (2011)]] for example. Some of these calculations are also available at [XrayDB Web App](#) (beta!).

### 1.3.1 X-ray attenuation by elements

The XrayDB database tabulates values of the X-ray mass attenuation coefficient,  $\mu/\rho$ , for each element. In most of the X-ray regime used in materials characterization (say, up to 150 keV), the photo-electric effect is the main process that causes X-ray attenuation. When the photo-electric process is dominant, the values for  $\mu/\rho$  depends strongly on  $Z$  of the atom and on X-ray energy  $E$ . In addition to these strong dependencies, sharp increases – so-called absorption edges – with be see at energies of bound core electron levels of atoms. To illustrate these characteristics, the following script will plot  $\mu/\rho$  for selected elements:

```
#!/usr/bin/env python
# XrayDB example script python/examples/mu_elements.py
#
# plot X-ray mass attenuation for selected elements
#
import numpy as np
import matplotlib.pyplot as plt
import wxmplot.interactive as wi
from xraydb import mu_elam , atomic_symbol

energy = np.arange(500, 120000, 10) # energy in eV

for elem in ('C', 'Cu', 'Au'):
    mu = mu_elam(elem, energy)
    plt.plot(energy, mu, label=elem, linewidth=2)

plt.title('X-ray mass attenuation')
plt.xlabel('Energy (eV)')
plt.ylabel(r'$\mu/\rho$ \rm\, (cm^2/gr)$')
plt.legend()
plt.yscale('log')
plt.xscale('log')
plt.show()
```

As you can see in Figure from this figure, the attenuation drops very strongly with  $E$  – approximately as  $E^3$ .  $\mu$  also depends strongly with  $Z$ , though the sharp absorption edges make this more complicated.

You can also observe that at relatively high energies for relatively low- $Z$  elements (such as C above about 20 keV) that the attenuation levels off. This is because the coherent (Rayleigh) and incoherent (Compton) scattering processes dominate, so that the photo-electric absorption is no longer the dominant X-ray scattering process. This can be illustrated by plotting the different components of  $\mu/\rho$  for C, as with the following script:

```
#!/usr/bin/env python
# XrayDB example script python/examples/mu_components_C.py
#
# plot components of X-ray mass attenuation for C
#
import numpy as np
```

(continues on next page)




Fig. 1: X-ray mass attenuation coefficient for C, Cu, and Au.

(continued from previous page)

```
import matplotlib.pyplot as plt
from xraydb import mu_elam

energy = np.arange(500, 120000, 10) # energy in eV

elem = 'C'
mu_total = mu_elam(elem, energy, kind='total')
mu_photo = mu_elam(elem, energy, kind='photo')
mu_incoh = mu_elam(elem, energy, kind='incoh')
mu_coher = mu_elam(elem, energy, kind='coh')

plt.title('X-ray mass attenuation for %s' % elem)
plt.plot(energy, mu_total, linewidth=2, label='Total')
plt.plot(energy, mu_photo, linewidth=2, label='Photo-electric')
plt.plot(energy, mu_incoh, linewidth=2, label='Incoherent')
plt.plot(energy, mu_coher, linewidth=2, label='Coherent')

plt.xlabel('Energy (eV)')
plt.ylabel(r'$\mu/\rho$ (cm2/gr)')
plt.legend()
plt.yscale('log')
plt.show()
```

which will generate the following plot:

Note that above 20 keV, the photo-electric absorption and incoherent Compton contributions are about equal, and that the Compton scattering dominates above 50 keV. As shown above, the photo-electric scattering will be much higher for heavier elements. The Rayleigh and Compton scattering have a much weaker dependence on  $Z$ , so that the photo-electric process dominates to higher energies. Replacing 'C' with 'Fe' in the script above will generate the following




Fig. 2: X-ray scattering and attenuation factors for C.

plot:

which shows that the Compton scattering reaching about 0.1 to 0.25 cm<sup>2</sup>/gr for Fe, about the same value as it was for C, while the photo-electric cross-section dominates past 100 keV.

### 1.3.2 $\mu$ calculations for materials

While one can use the above values for  $\mu/\rho$  to calculate the attenuation of X-rays by multi-element materials, the `material_mu()` function is available to do the more convenient calculation of the X-ray absorption coefficient  $\mu$  in units of 1/cm for a material and energy value and density (which are known for several common materials). This gives the length for which X-ray intensity is reduced by a factor of  $e$ , and so can be used to calculate the fraction of the X-rays transmitted through a material of known thickness, as  $\exp(-t\mu)$  for a material of thickness  $t$ . As a first example, we calculate the the fraction of X-ray transmitted through 1 mm of the water as a function of X-ray energy:

```
#!/usr/bin/env python
# XrayDB example script    python/examples/mu_water.py
#
# calculate the fraction of X-rays transmitted through 1 mm of water
#
import numpy as np
import matplotlib.pyplot as plt

from xraydb import material_mu

energy = np.linspace(1000, 41000, 201)

mu = material_mu('H2O', energy)

# mu is returned in 1/cm
```

(continues on next page)




Fig. 3: X-ray scattering and attenuation factors for Fe.

(continued from previous page)

```
trans = np.exp(-0.1 * mu)

plt.plot(energy, trans, label='transmitted')
plt.plot(energy, 1-trans, label='attenuated')
plt.title('X-ray absorption by 1 mm of water')
plt.xlabel('Energy (eV)')
plt.ylabel('Transmitted / Attenuated fraction')
plt.legend()
plt.show()
```

replacing:

```
mu = material_mu('H2O', energy)
```

with:

```
mu = material_mu('CaCO3', energy, density=2.71)
```

would generate the following plot

For many X-ray experiments, selecting the size of a material size so that its thickness is approximately 1 to 2 absorption length is convenient so that X-ray scattering and emission can be observed strongly, with neither all primary and scattered X-rays being absorbed by the material itself, but also not simply passing through the material without any interaction. For example, one can simply do:

```
>>> from xraydb import material_mu
>>> mu_20kev = xraydb.material_mu('CaCO3', 20000, density=2.71)
>>> print("CaCO3 1/e depth at 20keV = {:.3f} mm".format(10/mu_20kev))
CaCO3 1/e depth at 20keV = 0.648 mm
```




Fig. 4: Fraction of X-rays absorbed and transmitted by water




Fig. 5: Fraction of X-rays absorbed and transmitted by calcite

### 1.3.3 X-ray flux calculations for ionization chambers and photodiodes

Gas-filled ionization chambers are widely used as X-ray detectors. They are simple to use, inexpensive, and can be highly linear in estimating the photon flux over many orders of magnitude. X-rays entering a chamber filled with an inert gas (typically He, N<sub>2</sub>, or one of the noble gases, or a mixture of these) will be partially absorbed by the gas, with the strong energy dependence shown above. By adjusting the composition of the gas, nearly any fraction of the incident X-ray beam can be absorbed at a particular X-ray energy, making these ideal detectors to sample the intensity of an X-ray beam incident on a sample, while attenuating only a fraction of the beam.

Some of the X-rays in the gas will be absorbed by the photo-electric effect which will *ionize* the gas, generating free electrons and energetic ions. In fact, the first ionization will generate an electron-ion pair with the energy of the X-ray minus the binding energy of the core electron, but the high-energy electron and ion will further ionize other gas molecules. With an electric potential (typically on the order of 1 kV /cm) across the chamber, a current can be measured that is proportional to the X-ray energy and fluence of the X-rays.

In addition to the photo-electric absorption, X-rays can be attenuated by gas molecules in an ion chamber by incoherent (Compton) or coherent (Rayleigh) scattering processes. The coherent scattering will not generate any electrons in the gas, but will scatter X-rays out of the main beam. Incoherent scattering will generate some current, though only a small portion of the X-ray energy is given to an electron to generate a current. Compton scattering gives a distribution of energies to the scattered electron depending on the angle of scattering. The median energy (and a decent approximation of the mean energy) of electrons generated by Compton scattering X-rays of energy  $E$  will be

$$E_c = E / (1 + m_e c^2 / E)$$

For X-rays of 10 keV,  $E_c$  is about 192 eV. For 20 keV X-rays, it will be 750 eV, and for 50 keV X-rays, it will be 4.5 keV.

Although the energy is much less than the X-ray energy that would be fully converted to current in the photo-electric process the contribution can be important. This is especially true for low-Z gas molecules such as He and N<sub>2</sub> at relatively high energies (10 keV and above) for which incoherent scattering becomes much more important than photo-electric absorption, as shown above for C. This effect needs to be included for accurate estimates of fluxes from ion chamber currents.

The process of converting the X-ray generated current into X-ray fluence involves several steps. The energy from a single X-ray-generated electron is converted into a number of electron-ion pairs given by the *effective ionization potential* of the gas. These are reasonably well-known values (see [[Knoll (2010)]] that are all between 20 and 40 eV, given in the *Table of Effective Ionization Potentials*.

Table of Effective Ionization Potentials. Many of these are taken from [[Knoll (2010)]], while others appear to come from International Commission on Radiation Units & Measurement, Report 31, 1979. The names given are those supported by the functions `ionization_potential()` and `ionchamber_fluxes()`.

gas/material name(s)	potential (eV)
hydrogen, H	36.5
helium, He	41.3
nitrogen, N, N <sub>2</sub>	34.8
oxygen, O, O <sub>2</sub>	30.8
neon, Ne	35.4
argon, Ar	26.4
krypton, Kr	24.4
xenon, Xe	22.1
air	33.8
methane, CH <sub>4</sub>	27.3
carbondioxide, CO <sub>2</sub>	33.0
silicon, Si	3.68
germanium, Ge	2.97

From this table, we can see that the absorption (by photo-electric effect) of 1 X-ray of energy 10 keV will eventually generate about 300 electron-ion pairs. That is not much current, but if  $10^8$  Hz X-rays are absorbed per second, then the current generated will be around 5 nA. Of course, the thickness of the gas or more precisely the length of gas under ionizing potential will have an impact on how much current is generated. The photo-current will then be amplified and converted to a voltage using a current amplifier, and that voltage will then be recorded by a number of possible means. Note that while the ion chamber itself will be linear over many orders of magnitude of X-ray flux (provided the potential between the plates is high enough - typically in the 1 kV/cm range), a current amplifier at a particular setting of sensitivity will be linear only over a couple orders of magnitude (typically between output voltage of 0.05 to 5 V). Because of this, the sensitivity of the current amplifier used with an ion chamber needs careful attention.

A photo-diode works in much the same way as an ionization chamber. X-rays incident on the diode (typically Si or Ge) will be absorbed and generate a photo-current that can be collected. Typically PIN diodes are used, and with a small reverse bias voltage. Because the electrons do not need to escape the material but generate a current transported in the semiconductor, the effective ionization potential is much lower - a few times the semiconductor band gap instead of a few times the lowest core-level ionization potential. The current generated per X-ray will be larger than for an ion chamber, but still amplified with a current amplifier in the same way as is used for an ion chamber. Generally, diodes are thick enough that they absorb all incident X-rays.

The function `ionchamber_fluxes()` will help generate X-ray fluxes associated with an ion chamber and help handle all of these subtle issues, using the following inputs:

- *gas*: the gas, or mixture of gases used or 'Si' or 'Ge' for diodes.
- *length*: the length of the ion chamber, in cm.
- *energy*: the X-ray energy, in eV.
- *volts*: the output voltage of the current amplifier
- *sensitivity* and *sensitivity\_units*: the sensitivity or gain of the amplifier used to convert the photo-current to the recorded voltage.

The default *sensitivity\_units* is 'A/V' but can be set to any of the common SI prefixes such as 'p', 'pico', 'n', 'nano',  $\mu$ , (unicode '03bc'), 'u', 'micro', 'm', or 'milli', so that:

```
>>> fluxes = ionchamber_fluxes('N2', volts=1, energy=10000, length=10,
                                sensitivity=1.e-9)
>>> fluxes = ionchamber_fluxes('N2', volts=1, energy=10000, length=10,
                                sensitivity=1, sensitivity_units='nA/V')
```

will give the same results.

The output from `ionchamber_fluxes()` is a named tuple with 4 fields:

- *photo* - the flux absorbed by the photo-electric effect, in Hz.
- *incoherent* - the flux scattered by the Compton effects, in Hz.
- *incident* - the flux incident on the ion chamber, in Hz.
- *transmitted* - the flux beam leaving the ion chamber, in Hz.

Note that the ion chamber (on photo-diode) current is generated by the photo-electric and incoherent effects, which dominates for heavy elements and relatively low X-ray energies, but does not necessarily dominate at high X-ray energies. The photo-electric cross-section with the incident X-ray energy and the incoherent cross-section with the median Compton-scattering energy above are used to estimate the incident flux from the photo-current. The total attenuation cross-section, including the coherent cross-sections, is used to calculate the transmitted flu from the incident flux.

As an example calculation of ion chamber currents:

```
>>> fl = ionchamber_fluxes(gas='nitrogen', volts=1.25, energy=18000,
                           length=10.0, sensitivity=1.e-6)
>>> print(f"Incident= {fl.incident:g} Hz, Transmitted flux= {fl.transmitted:g} Hz")
Incident= 2.2358e+12 Hz, Transmitted flux= 2.214e+12 Hz
```

It is not uncommon for an ion chamber to be filled with a mixture of 2 or more gases so as to better control the fraction of X-rays absorbed in a chamber of fixed length. This can be specified by passing in a dictionary of gas name and fractional density, as with:

```
>>> fl = ionchamber_fluxes(gas={'Kr':0.5, 'Ar': 0.5}, volts=1.25,
                           energy=18000, length=10,
                           sensitivity=1, sensitivity_units='microA/V')
>>> print(f"Incident= {fl.incident:g} Hz, Transmitted flux= {fl.transmitted:g} Hz")
Incident= 1.43737e+10 Hz, Transmitted flux= 3.28986e+09 Hz
```

Finally, the pressure of the gas is sometimes adjusted to alter the fraction of the beam absorbed. The calculations here all use the densities at STP, but changes in gas density will be exactly linear to changing the length of the ion chamber.

### 1.3.4 X-ray mirror reflectivities

At very shallow angles of incidence X-rays can be reflected by total external reflection from a material. The reflectivity can be very high at relatively low energies and shallow angles, but drops off dramatically with increasing energy, increasing angle, and decreasing electron density. Still, this reflectivity is one of the few ways to steer X-ray beams and so is widely used in synchrotron radiation sources.

The reflectivity can be calculated with the `mirror_reflectivity()` function which takes X-ray energy, incident angle, and mirror material as arguments.

An example script, comparing the energy-dependence of the reflectivity for a few common mirror materials is given as

```
import numpy as np
from xraydb import mirror_reflectivity
import matplotlib.pyplot as plt

energy = np.linspace(1000, 51000, 501)

r_si = mirror_reflectivity('Si', 0.002, energy)
r_ni = mirror_reflectivity('Ni', 0.002, energy)
r_rh = mirror_reflectivity('Rh', 0.002, energy)
r_pt = mirror_reflectivity('Pt', 0.002, energy)

plt.plot(energy, r_si, label='Si')
plt.plot(energy, r_ni, label='Ni')
plt.plot(energy, r_rh, label='Rh')
plt.plot(energy, r_pt, label='Pt')

plt.title('X-ray reflectivity at  $\theta=2 \text{ mrad}$ ')
plt.xlabel('Energy (eV)')
plt.ylabel('Reflectivity')
plt.legend()
plt.show()
```






Fig. 6: X-ray mirror reflectivity at  $\theta = 2\text{mrad}$  for selected mirror surfaces and coatings used for mirrors.

### 1.3.5 Darwin widths of monochromator crystals

Bragg's law describes X-ray diffraction from crystals as

$$m\lambda = 2d \sin(\theta)$$

where  $\lambda$  is the X-ray wavelength,  $d$  the d-spacing of the crystal lattice plane,  $\theta$  the incident angle, and  $m$  the order of the reflection. For imperfect crystals, in which the lattice planes are not stacked perfectly over extended distances, the angular width of any particular reflection is dominated by the spread in d-spacing and the mosaicity inherent in the crystal. For perfect crystals, however, the angular width of a reflection is dominated by the fact that effectively all of the X-rays will scatter from the lattice well before any attenuation of the X-ray beam occurs. This *dynamical* diffraction gives a small but finite offset from the Bragg angle, and gives a broadened angular width to reflection. This is usually called the Darwin width (named for Charles G. Darwin, grandson of the more famous Charles R. Darwin). In addition, the refraction and in particular the absorption effects that give anomalous scattering (as calculated with `xray_delta_beta()`) make the “rocking curve” of reflected intensity as a function of angle an asymmetric shape.

All of these effects are included in the `darwin_width()` function, which follows very closely the description from chapter 6.4 in [[Als-Nielsen and McMorrow (2011)]]. The function takes inputs of

- *energy*: the X-ray energy, in eV.
- *crystal*: the atomic symbol for the crystal: ‘Si’, ‘Ge’, or ‘C’. [‘Si’]
- *hkl*: a tuple with (h, k, l) of the reflection used. [(1, 1, 1)]
- *a*: lattice constant [None - use nominal value for crystal]
- *polarization*: s, p, or u to specify the X-ray polarization relative to the crystal [s]
- *m*: the order of the reflection. [1]
- *ignore\_fl*: whether to ignore *fl*. [False]

- *ignore\_f2*: whether to ignore *f2*. [False]

Polarization of *s* should be used for vertically deflecting monochromators at most synchrotron sources (which will normally be horizontally polarized), and *p* should be used for horizontally deflecting monochromators. For crystals used to analyzed unpolarized X-ray emission, use *u*, which will give the average of *s* and *p* polarization.

As with *ionchamber\_fluxes()*, the output here is complicated enough that it is put into a named *DarwinWidth* tuple that will contain the following fields:

- *theta* - the nominal Bragg angle, in rad
- *theta\_offset* - the offset from the nominal Bragg angle, in rad.
- *theta\_width* - estimated angular Darwin width, in rad
- *theta\_fwhm* - estimated FWHM of the angular reflectivity curve, in rad
- *energy\_width* - estimated energy Darwin width, in eV
- *energy\_fwhm* - estimated FWHM energy reflectivity curve, in eV
- *zeta* - nd-array of  $\zeta = \Delta\lambda/\lambda$ .
- *dtheta* - nd-array of angles around from Bragg angle, in rad
- *denergy* - nd-array of energies around from Bragg energy, in eV
- *intensity* - nd-array of reflected intensity at *zeta* values.

Here, *dtheta* will be given by  $\Delta\theta = \zeta \tan(\theta)$ , and *denergy* will be given by  $\Delta E = \zeta E$ . All of the nd-arrays will be the same size, so that plots of reflectivity can be readily made. An example usage, printing the predicted energy and angular widths and plotting the intensity profile or “rocking curve” is

```
import numpy as np
from xraydb import darwin_width
import matplotlib.pyplot as plt

dw_si111 = darwin_width(10000, 'Si', (1, 1, 1))
dw_si333 = darwin_width(30000, 'Si', (3, 3, 3))

fmt_string = "Darwin Width for {:s} at {:.0f} keV: {:.5.2f} microrad, {:.5.2f} eV"
print(fmt_string.format('Si(111)', 10,
                        dw_si111.theta_width*1e6,
                        dw_si111.energy_width))

print(fmt_string.format('Si(333)', 30,
                        dw_si333.theta_width*1e6,
                        dw_si333.energy_width))

dtheta = dw_si111.dtheta*1e6
denergy = dw_si111.denergy[:-1]

# slightly advanced matplotlib hackery:
fig, ax = plt.subplots(constrained_layout=True)

ax.plot(dtheta, dw_si111.intensity, label='$I$, Si(111)', linewidth=2)
ax.plot(dtheta, dw_si111.intensity**2, label='$I^2$, Si(111)', linewidth=2)
ax.plot(dw_si333.dtheta*1e6, dw_si333.intensity**2, label='$I^2$ Si(333) 30 keV',
        linewidth=2)
```

(continues on next page)

(continued from previous page)

```

ax.set_title('X-ray diffraction intensity at 10keV')
ax.set_xlabel('Angle -  $\theta_B$  ( $\mu\text{rad}$ )')
ax.set_ylabel('Reflectivity')
ax.legend()

plt.show()

```

which will print out values of:

```

Darwin Width for Si(111) at 10 keV: 26.96 microrad, 1.34 eV
Darwin Width for Si(333) at 30 keV: 1.81 microrad, 0.27 eV

```

and generates a plot of




Fig. 7: X-ray monochromator diffracted intensities around the Si(111) reflection. Here,  $i$  represents the intensity of a single reflection, and  $i^2$  the intensity from 2 bounces, as for a double-crystal monochromator. The intensity and angular offset of the third harmonic is also shown.

Note that the values reported for *theta\_fwhm* and *energy\_fwhm* will be about 6% larger than the reported values for *theta\_width* and *energy\_width*. The *width* values closely follow the region of the curve where the reflectivity ignoring absorption would be 1 - the flat top of the curve. Since a double-crystal monochromator will suppress the tails of the reflectivity, this smaller value is the one typically reported as “the Darwin width”, though some sources will report this smaller value as “FWHM”.

## 1.4 Using XrayDB from Python

The *python* directory contains the source code for a Python module for XrayDB. This module gives a user-friendly wrapping of the XrayDB, and automates the the conversion of data from sqlite database into Python and numpy arrays. The module requires the *numpy*, *scipy* and *sqlalchemy* modules, all of which are readily available and can be installed with:

```
pip install xraydb
```

The current version of the Python module is 4.4.6, corresponding to version 6 of *xraydb.sqlite*.

### 1.4.1 The Python xraydb module

To use the XrayDB from Python, you can import the *xraydb* module and start using it:

```
>>> import xraydb
>>> xraydb.atomic_number('Ag')
47
#
# X-ray elastic (Thomson) scattering factors:
>>> import numpy as np
>>> q = np.linspace(0, 5, 11)
>>> xraydb.f0('Fe', q)
array([25.994603 , 11.50848469,  6.55945765,  4.71039413,  3.21048827,
        2.20939146,  1.65112769,  1.36705887,  1.21133507,  1.10155689,
        1.0035555 ])
```

```
#
# X-ray emission lines:
>>> for name, line in xraydb.xray_lines('Zn', 'K').items():
...     print(name, ' = ', line)
...
Ka3 = XrayLine(energy=8462.8, intensity=0.000316256, initial_level='K', final_level=
↳ 'L1')
Ka2 = XrayLine(energy=8614.1, intensity=0.294353, initial_level='K', final_level='L2
↳ ')
Ka1 = XrayLine(energy=8637.2, intensity=0.576058, initial_level='K', final_level='L3
↳ ')
Kb3 = XrayLine(energy=9567.6, intensity=0.0438347, initial_level='K', final_level=
↳ 'M2')
Kb1 = XrayLine(energy=9570.4, intensity=0.0846229, initial_level='K', final_level=
↳ 'M3')
Kb5 = XrayLine(energy=9648.8, intensity=0.000815698, initial_level='K', final_level=
↳ 'M4,5')
#
# X-ray absorption edges:
>>> xraydb.xray_edge('As', 'K')
XrayEdge(energy=11867.0, fyield=0.548989, jump_ratio=7.314)
#
# X-ray attenuation factors:
>>> as_kedge = xraydb.xray_edge('As', 'K').energy
>>> energies = np.linspace(-50, 50, 5) + as_kedge
>>> muvals = xraydb.mu_elam('As', energies)
>>> for en, mu in zip(energies, muvals):
...     print("{:.0f}    {:.2f}".format(en, mu))
...
11817    26.07
```

(continues on next page)

(continued from previous page)

11842	25.92
11867	25.77
11892	178.32
11917	177.38

### Table of XrayDB function for Atomic and X-ray data for the elements

Most of these function return some element-specific property from the element symbol or atomic number. Some of the data extends to Z=98 (Cf), but some data may not be available for Z > 92 (U). Except where noted, the data comes from [\[\[Elam, Ravel, and Sieber \(2002\)\]\]](#).

xraydb functions	description
<code>atomic_number()</code>	atomic number from symbol
<code>atomic_symbol()</code>	atomic symbol from number
<code>atomic_mass()</code>	atomic mass
<code>atomic_density()</code>	density of pure element
<code>f0()</code>	elastic scattering factor ( <a href="#">[[Waasmaier and Kirfel (1995)]]</a> )
<code>f0_ions()</code>	list of valid “ions” for <code>f0()</code> ( <a href="#">[[Waasmaier and Kirfel (1995)]]</a> )
<code>xray_edge()</code>	xray edge data for a particular element and edge
<code>xray_edges()</code>	dictionary of all X-ray edges data for an element
<code>xray_lines()</code>	dictionary of all X-ray emission line data for an element
<code>fluor_yield()</code>	fluorescent yield for an X-ray emission line
<code>ck_probability()</code>	Coster-Kronig transition probability between two atomic levels
<code>mu_elam()</code>	absorption cross-section, photo-electric or total for an element
<code>coherent_cross_section_elam()</code>	coherent scattering cross-section for an element
<code>incoherent_cross_section_elam()</code>	incoherent scattering cross-section for an element
<code>chantler_energies()</code>	energies of tabulation for Chantler data ( <a href="#">[[Chantler (2000)]]</a> )
<code>f1_chantler()</code>	$f'(E)$ anomalous scattering factor ( <a href="#">[[Chantler (2000)]]</a> )
<code>f2_chantler()</code>	$f''(E)$ anomalous scattering factor ( <a href="#">[[Chantler (2000)]]</a> )
<code>mu_chantler()</code>	absorption cross-section ( <a href="#">[[Chantler (2000)]]</a> )
<code>guess_edge()</code>	guess element and edge from energy of absorption edge
<code>chemparse()</code>	parse a chemical formula to atomic abundances
<code>validate_formula()</code>	test whether a chemical formula can be parsed.
<code>get_materials()</code>	get a dictionary of known materials {name:(formula, density)}
<code>get_material()</code>	get a (formula, density) tuple for a material in the materials database
<code>find_material()</code>	get a material instance for a material in the materials database
<code>add_material()</code>	add a material to local materials database
<code>material_mu()</code>	absorption cross-section for a material at X-ray energies
<code>material_mu_components()</code>	dictionary of elemental components of <i>mu</i> for material
<code>xray_delta_beta()</code>	anomalous index of refraction for material and energy
<code>darwin_width()</code>	Darwin widths for monochromator crystals
<code>mirror_reflectivity()</code>	X-ray reflectivities for mirror materials (thick slab limit)
<code>ionization_potential()</code>	effective ionization potential for a gas, as for ion chambers
<code>ionchamber_fluxes()</code>	calculate fluxes from ion chamber voltages, gases, and sensitivities

### `get_xraydb()`

return instance of the XrayDB

**Returns** XrayDB

### Example

```
>>> import xraydb
>>> xdb = xraydb.get_xraydb()
```

## 1.4.2 Atomic Properties

**atomic\_number** (*element*)

*z* for element name

**Parameters** **element** (*str*) – atomic symbol

**Returns** atomic number

**atomic\_symbol** (*z*)

atomic symbol for atomic number

**Parameters** **z** (*int*) – atomic number

**Returns** atomic symbol

**atomic\_mass** (*element*)

molar mass for an element

**Parameters** **element** (*int*, *str*) – atomic number, atomic symbol for element

**Returns** atomic mass, in AMU

**atomic\_density** (*element*)

density (gr/cm<sup>3</sup>) for common for of an element

**Parameters** **element** (*int*, *str*) – atomic number, atomic symbol for element

**Returns** density in gm/cm<sup>3</sup>

## 1.4.3 Elastic Scattering Factors

**f0** (*ion*, *q*)

elastic X-ray scattering factor, f0(q), for an ion.

**Parameters**

- **ion** (*int* or *str*) – atomic number, atomic symbol or ionic symbol of scatterer
- **q** (*float*, *ndarray*) – Q value(s) for scattering

**Returns** scattering factor for each Q value

### Notes

1. from D. Waasmaier and A. Kirfel, Acta Cryst. A51 p416 (1995) and International Tables for Crystallography, Vol. C.
2. *ion* can be of the form: 26, Fe, Fe2+. For a full list of ions use *f0\_ions()*
3. elements supported are from Z = 1 to 98 ('H' to 'Cf')
4.  $q = \sin(\theta) / \lambda$ , where  $\theta$ =incident angle,  $\lambda$ =X-ray wavelength

**f0\_ions** (*element=None*)

list ion names supported in the f0() calculation from Waasmaier and Kirfel.

**Parameters** *element* (*None*, *int*, *str*) – scatterer

**Returns** list of strings for matching ion names

### Notes

if element is None, all 211 ions are returned.

## 1.4.4 X-ray Edges

**xray\_edge** (*element*, *edge*, *energy\_only=False*)

get x-ray absorption edge data for an element: (energy(in eV), fluorescence yield, jump ratio)

### Parameters

- **element** (*int*, *str*) – atomic number, atomic symbol for element
- **edge** (*str*) – iupac symbol of X-ray edge
- **energy\_only** (*bool*) – whether to return only the energy [False]

**Returns** XrayEdge namedtuple containing (energy, fluorescence\_yield, edge\_jump) or float of energy

**xray\_edges** (*element*)

**get dictionary of x-ray absorption edges:** energy(in eV), fluorescence yield, and jump ratio for an element.

**Parameters** *element* (*int*, *str*) – atomic number, atomic symbol for element

**Returns** dictionary of XrayEdge named tuples.

### Notes

1. The dictionary will have keys of edge (iupac symbol) and values containing an XrayEdge namedtuple containing (energy, fluorescence\_yield, edge\_jump)

**core\_width** (*element*, *edge=None*)

returns core hole width for an element and edge

### Parameters

- **element** (*int* or *str*) – element
- **edge** (*None* or *str*) – edge to consider

**Returns** core width or list of core widths

## Notes

1. if edge is None, values are return for all edges
2. Data from Krause and Oliver (1979) and Keski-Rahkonen and Krause (1974)

**guess\_edge** (*energy*, *edges*=['K', 'L3', 'L2', 'L1', 'M5'])  
guess an element and edge based on energy (in eV)

### Parameters

- **energy** (*float*) – approximate edge energy (in eV)
- **edges** (*None or list of strings*) – edges to consider

**Returns** a tuple of (atomic symbol, edge) for best guess

## Notes

by default, the list of edges is ['K', 'L3', 'L2', 'L1', 'M5']

## 1.4.5 X-ray Emission Lines

**xray\_lines** (*element*, *initial\_level*=None, *excitation\_energy*=None)  
get dictionary of X-ray emission lines of an element

### Parameters

- **element** (*int, str*) – atomic number, atomic symbol for element
- **initial\_level** (*None or str*) – iupac symbol of initial level
- **excitation\_energy** (*None or float*) – exciation energy

**Returns** dict of X-ray lines with keys of siegbahn notation and values of XrayLine tuples of (energy, intensity, initial level, final level)

## Notes

1. excitation energy will supercede initial\_level, as it means 'all intial levels with below this energy

### Example:

```
>>> for name, line in xraydb.xray_lines('Mn', 'K').items():  
...     print(name, line)  
...  
Ka3 XrayLine(energy=5769.9, intensity=0.000265963, initial_level='K', final_  
↪level='L1')  
Ka2 XrayLine(energy=5889.1, intensity=0.293941, initial_level='K', final_  
↪level='L2')  
Ka1 XrayLine(energy=5900.3, intensity=0.58134, initial_level='K', final_level=  
↪'L3')  
Kb3 XrayLine(energy=6491.8, intensity=0.042234, initial_level='K', final_  
↪level='M2')  
Kb1 XrayLine(energy=6491.8, intensity=0.0815329, initial_level='K', final_  
↪level='M3')  
Kb5 XrayLine(energy=6537.0, intensity=0.000685981, initial_level='K', final_  
↪level='M4,5')
```



**fluor\_yield** (*element, edge, line, energy*)

fluorescence yield for an X-ray emission line or family of lines.

**Parameters**

- **element** (*int, str*) – atomic number, atomic symbol for element
- **edge** (*str*) – iupac symbol of X-ray edge
- **line** (*str*) – siegbahn notation for emission line
- **energy** (*float*) – incident X-ray energy

**Returns** fluorescence yield, weighted average fluorescence energy, net\_probability

**Examples**

```
>>> xraydb.fluor_yield('Fe', 'K', 'Ka', 8000)
0.350985, 6400.752419799043, 0.874576096
```

```
>>> xraydb.fluor_yield('Fe', 'K', 'Ka', 6800)
0.0, 6400.752419799043, 0.874576096
```

```
>>> xraydb.fluor_yield('Ag', 'L3', 'La', 6000)
0.052, 2982.129655446868, 0.8618990000000000
```

**See also:**

*xray\_lines* which gives the full set of emission lines ('Ka1', 'Kb3', etc) and probabilities for each of these.

**ck\_probability** (*element, initial, final, total=True*)

transition probability for an element, initial, and final levels.

**Parameters**

- **element** (*int, str*) – atomic number, atomic symbol for element
- **initial** (*str*) – iupac symbol for initial level
- **final** (*str*) – iupac symbol for final level
- **total** (*bool*) – whether to include transitions via possible intermediate levels [True]

**Returns** transition probability, or 0 if transition is not allowed.

## 1.4.6 Absorption and Scattering Cross-sections

**mu\_elam** (*element, energy, kind='total'*)

X-ray mass attenuation coefficient, mu/rho, for an element and energy or array of energies. Data is from the Elam tables.

**Parameters**

- **element** (*int, str*) – atomic number, atomic symbol for element
- **energy** (*float or ndarray*) – energy or array of energies
- **kind** (*str*) – type of cross-section to use, one of ('total', 'photo', 'coh', 'incoh') ['total']

**Returns** float value or ndarray

### Notes

1. Values returned are in units of  $\text{cm}^2/\text{gr}$
2. The default is to return total attenuation coefficient.

**coherent\_cross\_section\_elam** (*element, energy*)

coherent scattering cross-section for an element and energy or array of energies. Data is from the Elam tables.

#### Parameters

- **element** (*int, str*) – atomic number, atomic symbol for element
- **energy** (*float or ndarray*) – energy or array of energies

**Returns** float value or ndarray

### Notes

1. Values returned are in units of  $\text{cm}^2/\text{gr}$

**incoherent\_cross\_section\_elam** (*element, energy*)

incoherent scattering cross-section for an element and energy or array of energies. Data is from the Elam tables.

#### Parameters

- **element** (*int, str*) – atomic number, atomic symbol for element
- **energy** (*float or ndarray*) – energy or array of energies

**Returns** float value or ndarray

### Notes

1. Values returned are in units of  $\text{cm}^2/\text{gr}$

**chantler\_energies** (*element, emin=0, emax=1000000000.0*)

energies at which Chantler data is tabulated for a particular element.

#### Parameters

- **element** (*int, str*) – atomic number, atomic symbol for element
- **emin** (*float*) – lower bound of energies (default=0)
- **emax** (*float*) – upper bound of energies (default=1.e9)

**Returns** ndarray of energies

## Notes

energies are in eV

**f1\_chantler** (*element, energy, \*\*kws*)

real part of anomalous x-ray scattering factor for an element and energy or array of energies. Data is from the Chantler tables.

### Parameters

- **element** (*int, str*) – atomic number, atomic symbol for element
- **energy** (*float or ndarray*) – energy or array of energies

**Returns** float value or ndarray

## Notes

1. Values returned are in units of electrons

**f2\_chantler** (*element, energy*)

imaginary part of anomalous x-ray scattering factor for an element and energy or array of energies. Data is from the Chantler tables.

### Parameters

- **element** (*int, str*) – atomic number, atomic symbol for element
- **energy** (*float or ndarray*) – energy or array of energies

**Returns** float value or ndarray

## Notes

1. Values returned are in units of electrons

**mu\_chantler** (*element, energy, incoh=False, photo=False*)

X-ray mass attenuation coefficient,  $\mu/\rho$ , for an element and energy or array of energies. Data is from the Chantler tables.

### Parameters

- **element** (*int, str*) – atomic number, atomic symbol for element
- **energy** (*float or ndarray*) – energy or array of energies
- **incoh** (*bool*) – whether to return only the incoherent contribution [False]
- **photo** (*bool*) – whether to return only the photo-electric contribution [False]

**Returns** float value or ndarray

## Notes

1. Values returned are in units of cm<sup>2</sup>/gr
2. The default is to return total attenuation coefficient.

## 1.4.7 Chemical and Materials database

### **chemparse** (*formula*)

parse a chemical formula to a dictionary of elemental abundances

**Parameters** **formula** (*str*) – chemical formula

**Returns** dict of element symbol and abundance.

### Examples

```
>>> from xraydb import chemparse
>>> chemparse('Mn(SO4)2(H2O)7')
{'H': 14.0, 'S': 2.0, 'Mn': 1, 'O': 15.0}
```

```
>>> chemparse('Zn1.e-5Fe3O4')
{'Zn': 1e-05, 'Fe': 3.0, 'O': 4.0}
```

```
>>> chemparse('CO')
{'C': 1, 'O': 1}
>>> chemparse('Co')
{'Co': 1}
```

```
>>> chemparse('co')
ValueError: unrecognized element or number:
co
```

### **validate\_formula** (*formula*)

return whether a chemical formula is valid and can be parsed to a dictionary with chemparse()

**Parameters** **formula** (*str*) – chemical formula

**Returns** bool (True or False) for whether chemparse() will succeed

### Examples

```
>>> from xraydb import validate_formula
>>> validate_formula('Mn(SO4)2(H2O)7')
True
```

```
>>> validate_formula('Mn(SO4)2(H2O7)')
False
```

```
>>> validate_formula('Z')
False
```

### **get\_materials** (*force\_read=False, categories=None*)

get dictionary of all available materials

**Parameters**

- **force\_read** (*bool*) – whether to force a re-reading of the materials database [False]
- **categories** (*list of strings or None*) – restrict results to those that match category names

**Returns** dict with keys of material name and values of Materials instances

**Examples**

```
>>> for name, m in xraydb.get_materials().items():
...     print(name, m)
...
water H2O 1.0
lead Pb 11.34
aluminum Al 2.7
kapton C22H10N2O5 1.42
polyimide C22H10N2O5 1.42
nitrogen N 0.00125
argon Ar 0.001784
...
```

**find\_material** (*name*)

look up material name, return material instance

**Parameters** **name** (*str*) – name of material or chemical formula

**Returns** material instance

**Examples**

```
>>> xraydb.find_material('kapton')
Material(formula='C22H10N2O5', density=1.42, name='kapton', categories=['polymer
↪'])
```

**See also:**

get\_material()

**get\_material** (*name*)

look up material name, return formula and density

**Parameters** **name** (*str*) – name of material or chemical formula

**Returns** chemical formula, density of material

**Examples**

```
>>> xraydb.get_material('kapton')
('C22H10N2O5', 1.43)
```

**See also:**

find\_material()

**add\_material** (*name, formula, density, categories=None*)

add a material to the users local material database

**Parameters**

- **name** (*str*) – name of material
- **formula** (*str*) – chemical formula
- **density** (*float*) – density
- **categories** (*list of strings or None*) – list of category names

**Returns** None**Notes**

the data will be saved to \$HOME/.config/xraydb/materials.dat in the users home directory, and will be useful in subsequent sessions.

**Examples**

```
>>> xraydb.add_material('becopper', 'Cu0.98e0.02', 8.3, categories=['metal'])
```

## 1.4.8 X-ray properties of materials

For some further examples, see *Example Calculations of X-ray properties of materials*.

**material\_mu** (*name, energy, density=None, kind='total'*)

X-ray attenuation length (in 1/cm) for a material by name or formula

**Parameters**

- **name** (*str*) – chemical formul or name of material from materials list.
- **energy** (*float or ndarray*) – energy or array of energies in eV
- **density** (*None or float*) – material density (gr/cm<sup>3</sup>).
- **kind** (*str*) – ‘photo’ or ‘total’ for whether to return the photo-absorption or total cross-section [‘total’]

**Returns** absorption length in 1/cm**Notes**

1. material names are not case sensitive, chemical compounds are case sensitive.
2. mu\_elam() is used for mu calculation.
3. if density is None and material is known, that density will be used.

## Examples

```
>>> material_mu('H2O', 10000.0)
5.32986401658495
```

**material\_mu\_components** (*name, energy, density=None, kind='total'*)  
 material\_mu\_components: absorption coefficient (in 1/cm) for a compound

### Parameters

- **name** (*str*) – chemical formul or name of material from materials list.
- **energy** (*float or ndarray*) – energy or array of energies in eV
- **density** (*None or float*) – material density (gr/cm<sup>3</sup>).
- **kind** (*str*) – ‘photo’ or ‘total’ for whether to return photo-absorption or total cross-section [‘total’]

### Returns

**dict for constructing mu per element**, with elements ‘mass’ (total mass), ‘density’, and ‘elements’ (list of atomic symbols for elements in material). For each element, there will be an item (atomic symbol as key) with tuple of (stoichiometric fraction, atomic mass, mu)

## Examples

```
>>> xraydb.material_mu('quartz', 10000)
50.36774553547068
>>> xraydb.material_mu_components('quartz', 10000)
{'mass': 60.0843, 'density': 2.65, 'elements': ['Si', 'O'],
 'Si': (1, 28.0855, 33.87943243018506), 'O': (2.0, 15.9994, 5.952824815297084)}
```

**xray\_delta\_beta** (*material, density, energy*)  
 anomalous components of the index of refraction for a material, using the tabulated scattering components from Chantler.

### Parameters

- **material** – chemical formula (‘Fe2O3’, ‘CaMg(CO3)2’, ‘La1.9Sr0.1CuO4’)
- **density** – material density in g/cm<sup>3</sup>
- **energy** – x-ray energy in eV

**Returns** (delta, beta, atlen)

**where** delta : real part of index of refraction beta : imag part of index of refraction atlen : attenuation length in cm

These are the anomalous scattering components of the index of refraction:

$$n = 1 - \delta - i \cdot \beta = 1 - \frac{\lambda^2}{2\pi} \sum_j (n_j \cdot f_j)$$

Adapted from code by Yong Choi

**darwin\_width** (*energy, crystal='Si', hkl=(1, 1, 1), a=None, polarization='s', ignore\_f2=False, ignore\_f1=False, m=1*)  
 darwin width for a crystal reflection and energy

### Parameters

- **energy** (*float*) – X-ray energy in eV
- **crystal** (*string*) – name of crystal (one of ‘Si’, ‘Ge’, or ‘C’) [‘Si’]
- **hkl** (*tuple*) – h, k, l for reflection [(1, 1, 1)]
- **a** (*float or None*) – lattice constant [None - use built-in value]
- **polarization** (*'s', 'p', 'u'*) – mono orientation relative to X-ray polarization [‘s’]
- **ignore\_f1** (*bool*) – ignore contribution from f1 - dispersion (False)
- **ignore\_f2** (*bool*) – ignore contribution from f2 - absorption (False)
- **m** (*int*) – order of reflection [1]

### Returns

A named tuple ‘DarwinWidth’ with the following fields

*theta*: float, nominal Bragg angle, in rad,  
*theta\_offset*: float, angular offset from Bragg angle, in rad,  
*theta\_width*: float, estimated angular Darwin width, in rad,  
*theta\_fwhm*: float, estimated FWHM of angular intensity, in rad,  
*energy\_width*: float, estimated angular Darwin width, in rad,  
*energy\_fwhm*: float, estimated FWHM of energy intensity, in eV,  
*zeta*: nd-array of Zeta parameter ( $\Delta\lambda / \lambda$ ),  
*dtheta*: nd-array of angles away from Bragg angle, theta in rad,  
*denergy*: nd-array of energies away from Bragg energy, in eV,  
*intensity*: nd-array of reflected intensity

### Notes

1. This follows the calculation from section 6.4 of Elements of Modern X-ray Physics, 2nd Edition J Als-Nielsen, and D. McMorrow.
2. Only diamond structures (Si, Ge, diamond) are currently supported. Default values of lattice constant *a* are in Angstroms: 5.4309 for Si, 5.6578, for ‘Ge’, and 3.567 for ‘C’.
3. The *theta\_width* and *energy\_width* values will closely match the width of the intensity profile that would = 1 when ignoring the effect of absorption. These are the values commonly reported as ‘Darwin Width’. The value reported for *theta\_fwhm* and *energy\_fwhm* are larger than this by  $\sqrt{9/8} \approx 1.06$ .
4. Polarization can be ‘s’, ‘p’, ‘u’, or None. ‘s’ means vertically deflecting crystal and a horizontally-polarized source, as for most synchrotron beamlines. ‘p’ is for a horizontally-deflecting crystal. ‘u’ or None is for unpolarized light, as for most fluorescence/emission.



## Examples

```
>>> dw = darwin_width(10000, crystal='Si', hkl=(1, 1, 1))
>>> dw.theta_width, dw.energy_width
(2.8593683930207114e-05, 1.4177346002236872)
```

**mirror\_reflectivity** (*formula, theta, energy, density=None, roughness=0.0, polarization='s'*)  
mirror reflectivity for a thick, singl-layer mirror.

### Parameters

- **formula** (*string*) – material name or formula ('Si', 'Rh', 'silicon')
- **theta** (*float or nd-array*) – mirror angle in radians
- **energy** (*float or nd-array*) – X-ray energy in eV
- **density** (*float or None*) – material density in g/cm<sup>3</sup>
- **roughness** (*float*) – mirror roughness in Angstroms
- **polarization** (*'s' or 'p'*) – mirror orientation relative to X-ray polarization

**Returns** mirror reflectivity values

### Notes

1. only one of theta or energy can be an nd-array
2. density can be *None* for known materials
3. polarization of 's' puts the X-ray polarization along the mirror surface, 'p' puts it normal to the mirror surface. For horizontally polarized X-ray beams from storage rings, 's' will usually mean 'vertically deflecting' and 'p' will usually mean 'horizontally deflecting'.

**ionization\_potential** (*gas*)

return effective ionization potential for a gas or diode semiconductor, as appropriate for ionization chambers in the linear regime (not in the 'proportional counter' regime) or for PIN photodiodes (not in 'avalanche' mode).

**Parameters** **gas** (*string*) – name of gas or 'Si' or 'Ge'

**Returns** ionization potential in eV

### Notes

Data from G. F. Knoll, Radiation Detection and Measurement, Table 5-1, and from ICRU Report 31, 1979. Supported gas names and effective potentials:

gas names	potential (eV)
hydrogen, H	36.5
helium, He	41.3
nitrogen, N, N2	34.8
oxygen, O, O2	30.8
neon, Ne	35.4
argon, Ar	26.4
krypton, Kr	24.4
xenon, Xe	22.1
air	33.8
methane, CH4	27.3
carbondioxide, CO2	33.0
silicon, Si	3.68
germanium, Ge	2.97

If the gas is not recognized the default value of 32.0 eV will be returned.

**ionchamber\_fluxes** (*gas*='nitrogen', *volts*=1.0, *length*=100.0, *energy*=10000.0, *sensitivity*=1e-06, *sensitivity\_units*='A/V', *with\_compton*=1)

return ion chamber and PIN diode fluxes for a gas, mixture of gases, or semiconductor material, ion chamber length (or diode thickness), X-ray energy, recorded voltage and current amplifier sensitivity. See note for details.

#### Parameters

- **gas** (*string* or *dict*) – name or formula of fill gas (see note 1) ['nitrogen']
- **volts** (*float*) – measured voltage output of current amplifier [1.0]
- **length** (*float*) – active length of ion chamber in cm [100]
- **energy** (*float*) – X-ray energy in eV [10000]
- **sensitivity** (*float*) – current amplifier sensitivity [1.e-6]
- **sensitivity\_units** (*string*) – units of current amplifier sensitivity (see note 2 for options) ['A/V']
- **with\_compton** (*int*) – switch to control the contribution of Compton scattering (see note 3) [1]

#### Returns

named tuple IonchamberFluxes with fields

*incident* flux of beam incident on ion chamber in Hz

*transmitted* flux of beam output of ion chamber in Hz

*photo* flux absorbed by photo-electric effect in Hz

*incoherent* flux attenuated by incoherent scattering in Hz

## Examples

```
>>> from xraydb import ionchamber_fluxes
>>> fl = ionchamber_fluxes(gas='helium', volts=1.25, length=20.0,
                           energy=10000.0, sensitivity=1.e-9)
```

```
>>> print(f"Fluxes: In={fl.incident:g}, Out={fl.transmitted:g}, Transmitted=
↳ {100*fl.transmitted/fl.incident:.2f}%")
Fluxes: In=2.79631e+11, Out=2.79383e+11, Transmitted=99.91%
```

```
>>> fl = ionchamber_fluxes(gas='nitrogen', volts=1.25, length=20.0,
                           energy=10000.0, sensitivity=1.e-6)
```

```
>>> print(f"Fluxes: In={fl.incident:g}, Out={fl.transmitted:g}, Transmitted=
↳ {100*fl.transmitted/fl.incident:.2f}%")
Fluxes: In=3.20045e+11, Out=2.90464e+11, Transmitted=90.76%
```

```
>>> fl = ionchamber_fluxes(gas={'nitrogen':0.5, 'helium': 0.5},
                           volts=1.25, length=20.0, energy=10000.0,
                           sensitivity=1.e-6)
```

```
>>> print(f"Fluxes: In={fl.incident:g}, Out={fl.transmitted:g}, Transmitted=
↳ {100*fl.transmitted/fl.incident:.2f}%")
Fluxes: In=6.83845e+11, Out=6.51188e+11, Transmitted=95.22%
```

## Notes

1. The gas value can either be a string for the name of chemical formula for the gas or diode material, or dictionary with keys that are gas names or formulas and values that are the relative fraction for mixes gases. For diode materials, mixtures are not supported.

The gas formula is used both the contributions for mu and to get the weighted effective ionization potential for the material.

The effective ionization potentials are known for a handful of gases and diodes (see *ionization\_potential* function), and range between 20 and 45 eV for gases, and around 3 eV for semiconductors. For unknown gases the value of 32.0 eV will be used.

2. The *sensitivity* and *sensitivity\_units* arguments have some overlap to specify the sensitivity or gain of the current amplifier. Generally, the units in A/V, but you can add a common SI prefix of 'p', 'pico', 'n', 'nano', (unicode 'u03bc'), 'u', 'micro', 'm', 'milli' so that, *ionchamber\_fluxes(..., sensitivity=1.e-6)* and *ionchamber\_fluxes(..., sensitivity=1, sensitivity\_units='uA/V')* will both give a sensitivity of 1 microAmp / Volt .
3. The effect of Compton scattering on the ion chamber current can be approximated in 3 ways with different values for *with\_compton* (default=1):
  - >0: use Compton-shifted electron energy (best approximation)
  - =0: completely ignore the effect of Compton scattering on current
  - <0: use incident energy as Compton energy.

Using -1 should reproduce the calculation from Hephaestus reasonably well, but using 1 is highly recommended.

4. Ion Chamber Current calculation: the total attenuation in an ion chamber includes photo, incoherent (Compton), and coherent (Rayleigh) scattering contributions:

$$\text{Flux\_transmitted} = \text{Flux\_in} * \exp(-t * \mu_{\text{total}})$$

The current in an Ion Chamber has a contribution from both the photo-electric and incoherent (Compton) cross-section, but not the coherent scattering. For an incident flux (in Hz) of  $\text{Flux\_in}$ , the fluxes of X-rays attenuated by the different processes are:

$$\text{Flux\_photo} = \text{Flux\_in} * [1 - \exp(-t * \mu_{\text{photo}})]$$

$$\text{Flux\_incoh} = \text{Flux\_in} * [1 - \exp(-t * \mu_{\text{incoh}})]$$

$$\text{Flux\_coh} = \text{Flux\_in} * [1 - \exp(-t * \mu_{\text{coh}})]$$

$$\text{Flux\_transmitted} = \text{Flux\_in} - \text{Flux\_total}$$

For  $\text{Flux\_photo}$ , all of the X-ray energy is converted to current:

$$\text{IC\_current\_photo} = \text{Flux\_photo} * \text{Energy} * q_e / \text{ion\_pot}$$

where  $q_e$  is the electron charge ( $1.602 \times 10^{-19}$  C), and  $\text{ion\_pot}$  is the effective ionization potential (see *ionization\_potential*).

For  $\text{Flux\_incoh}$ , the energy transferred to the electron is converted to electron and ion current as:

$$\text{IC\_current\_incoh} = \text{Flux\_incoh} * \text{Energy\_Compton} * q_e / \text{ion\_pot}$$

where:

$$\text{Energy\_Compton} = \text{Energy} / (1 + m_e c^2 / \text{Energy})$$

is the approximate energy of the Compton-scattered electron. This estimate of the energy is a simplification: the transferred energy will be angle-dependent, ranging from 0 to  $2 * \text{Energy\_Compton}$ , with a distribution that depends on Energy and polarization.

## 1.5 Overview of Atomic and X-ray Data

The data provided in XrayDB includes Atomic data and characteristic energies and cross sections for the interaction of X-rays with elements. A few definitions and conventions necessary for using this data are discussed here.

### 1.5.1 Elements

Most of the data resources are accessed by an elements *Atomic Symbol*. For the Python module, most methods will take *element* as the first argument, and this can either be the integer atomic number or the string for the atomic symbol.

### 1.5.2 Physical Units

Elemental densities are given in  $\text{gr/cm}^3$ , and molar masses are given in AMU. Unless otherwise stated, all energies are in units of eV.

### 1.5.3 X-ray Edges

Several resources (database tables, python methods) take either an `edge` or a `level` argument to signify a core electronic level. These are strings and must be one of the levels listed in the [Table of X-ray edge names](#).

Table of X-ray Edges and Core electronic levels. The Names are the IUPAC symbols for the core electronic levels.

Name	electronic level	Name	electronic level
K	1s	N5	4d <sub>5/2</sub>
L3	2p <sub>3/2</sub>	N4	4d <sub>3/2</sub>
L2	2p <sub>1/2</sub>	N3	4p <sub>3/2</sub>
L1	2s	N2	4p <sub>1/2</sub>
M5	3d <sub>5/2</sub>	N1	4s
M4	3d <sub>3/2</sub>	O3	5p <sub>3/2</sub>
M3	3p <sub>3/2</sub>	O2	5p <sub>1/2</sub>
M2	3p <sub>1/2</sub>	O1	5s
M1	3s	P3	6p <sub>3/2</sub>
N7	4f <sub>7/2</sub>	P2	6p <sub>1/2</sub>
N6	4f <sub>5/2</sub>	P1	6s

### 1.5.4 X-ray Lines

Many resources (database tables or methods) take emission line arguments. These are all strings and follow the latinized version of the Siegbahn notation as indicated in the [Table of X-ray emission line names](#).

Table of X-ray emission line names and the corresponding Siegbahn and IUPAC notations

Name	IUPAC	Siegbahn	Name	IUPAC	Siegbahn
Ka1	K-L3	$K\alpha_1$	Lb4	L1-M2	$L\beta_4$
Ka2	K-L2	$K\alpha_2$	Lb5	L3-O4,5	$L\beta_5$
Ka3	K-L1	$K\alpha_3$	Lb6	L3-N1	$L\beta_6$
Kb1	K-M3	$K\beta_1$	Lg1	L2-N4	$L\gamma_1$
Kb2	K-N2,3	$K\beta_2$	Lg2	L1-N2	$L\gamma_2$
Kb3	K-M2	$K\beta_3$	Lg3	L1-N3	$L\gamma_3$
Kb4	K-N4,5	$K\beta_2$	Lg6	L2-O4	$L\gamma_6$
Kb5	K-M4,5	$K\beta_3$	Ll	L3-M1	$Ll$
La1	L3-M5	$L\alpha_1$	Ln	L2-M1	$L\nu$
La2	L3-M4	$L\alpha_1$	Ma	M5-N6,7	$M\alpha$
Lb1	L2-M4	$L\beta_1$	Mb	M4-N6	$M\beta$
Lb2,15	L3-N4,5	$L\beta_2, L\beta_{15}$	Mg	M3-N5	$M\gamma$
Lb3	L1-M3	$L\beta_3$	Mz	M4,5-N6,7	$M\zeta$

## 1.5.5 Cross Sections

The photo-absorption and scattering cross sections from [\[Elam, Ravel, and Sieber \(2002\)\]](#) and [\[Chantler \(2000\)\]](#) are in  $\text{cm}^2/\text{gr}$ .

The data from [\[Elam, Ravel, and Sieber \(2002\)\]](#) is held as logarithms of energy, cross section, and logarithm of the 2nd derivative of cross section that allows for cubic spline interpolation in log-log space.

## 1.6 Using the XrayDB xraydb.sqlite

All the data for the X-ray database is held in the SQLite3 file `xraydb.sqlite`. To use with SQLite, this file is all you need. While many programs and languages can access SQLite files, basic usage with the `sqlite3` program (available from Windows, Mac OS X, and Linux) can be as simple as:

```
system~> sqlite3 xraydb.sqlite
sqlite> .headers on
sqlite> select * from elements where atomic_number=47;
atomic_number|element|molar_mass|density
47|Ag|107.868|10.48
```

That is, you can retrieve the data using standard SQL queries built-in to SQLite. Of course, the expectation is that you'd want to use this database within a programming environment. Currently, wrappers exist only for Python.

### 1.6.1 Overall Database Schema

The schema for the SQLite3 database describes where data is held in the database, and how to access it. The schema for the current version (4) looks like this:

```
Table Version (id integer primary key,
               tag text,
               date text,
               notes text);
Table elements (atomic_number integer primary key,
                element text,
                molar_mass real,
                density real);
Table xray_levels (id integer primary key,
                  element text,
                  iupac_symbol text,
                  absorption_edge real,
                  fluorescence_yield real,
                  jump_ratio real);
Table xray_transitions (id integer primary key,
                       element text,
                       iupac_symbol text,
                       siegbahn_symbol text,
                       initial_level text,
                       final_level text,
                       emission_energy real,
                       intensity real);
Table Coster_Kronig (id integer primary key,
                    element text,
                    initial_level text,
                    final_level text,
```

(continues on next page)

(continued from previous page)

```

        transition_probability real,
        total_transition_probability real);
Table photoabsorption (id integer primary key,
        element text,
        log_energy text,
        log_photoabsorption text,
        log_photoabsorption_spline text);
Table scattering (id integer primary key,
        element text,
        log_energy text,
        log_coherent_scatter text,
        log_coherent_scatter_spline text,
        log_incoherent_scatter text,
        log_incoherent_scatter_spline text);
Table Waasmaier (id integer primary key,
        atomic_number integer,
        element text,
        ion text,
        offset real,
        scale text,
        exponents text);
Table KeskiRahkonen_Krause (id integer primary key,
        atomic_number integer,
        element text,
        edge text,
        width float);
Table Krause_Oliver (id integer primary key,
        atomic_number integer,
        element text,
        edge text,
        width float);
Table corelevel_widths (id integer primary key,
        atomic_number integer,
        element text,
        edge text,
        width float);
Table Chantler (id integer primary key,
        element text,
        sigma_mu real,
        mue_f2 real,
        density real,
        corr_henke float,
        corr_cl35 float,
        corr_nucl float,
        energy text,
        f1 text,
        f2 text,
        mu_photo text,
        mu_incoh text,
        mu_total text);

```

More details for each table are given below.

**Note:** in the tables below the type of *json array* means that arrays of numerical data are stored in the database as text of JSON-encoded arrays.

### 1.6.2 Version Table

The *Version* table holds data about the revisions to the database file itself. Each row represents a single revision.

DB Table of Database Versions

Column	Type	Description
id	integer	counter (primary tag)
tag	text	version name
date	text	date string
notes	text	notes on changes for version

### 1.6.3 Elements Table

The *elements* table holds basic data about each element. Each row represents an element.

DB Table of Basic Properties of the Elements

Column	Type	Description
atomic_number	integer	Atomic Number, Z
element	text	Atomic symbol
molar_mass	float	Atomic mass in AMU
density	float	Density of pure element (gr/cm <sup>3</sup> )

### 1.6.4 Xray\_Levels Table

The *xray\_levels* table holds data for electronic levels of atoms. Each row represents a core electronic level.

DB Table of X-ray and core electronic levels. *fluorescence yield* gives the probability of an empty level refilling by X-ray fluorescence. The *jump ratio* is the ratio of values for photo-electric cross section (that is, from [Photoabsorption Table](#)) 1 eV above the absorption edge to that 1 eV below the absorption edge. See [Table of X-ray Edges](#)

Column	Type	Description
id	integer	Index (primary key)
element	text	Atomic symbol for element
iupac_symbol	text	IUPAC symbol for level ('K','L3',...)
absorption_edge	float	binding energy for level (eV)
fluorescence_yield	float	fluorescence yield (fraction)
jump_ratio	float	ratio of mu_photo across edge

### 1.6.5 Xray\_Transitions Table

The *xray\_transitions* table holds data for transitions between electronic levels of atoms. Each row represents a transition between two levels.

DB Table of X-ray Transitions. Both IUPAC and Siegbahn symbols are given (see [Table of X-ray emission lines](#)), as well as the initial and final levels. The *intensity* is the relative intensity of the transition for a given *initial level*.



Column	Type	Description
id	integer	Index (primary key)
element	text	Atomic symbol for element
iupac_symbol	text	IUPAC symbol for transition
siegbahn_symbol	text	Siegbahn symbol for transition
initial_level	text	IUPAC symbol for initial level
final_level	text	IUPAC symbol for final level
emission_energy	float	fluorescence energy (eV)
intensity	float	relative intensity for transition

### 1.6.6 Photoabsorption Table

The *photoabsorption* table holds data for the photo-electric absorption cross sections in  $\text{cm}^2/\text{gr}$ . Each row represents an element.

DB Table of Photoabsorption Cross Sections. JSON-encoded arrays are held for logs of energy, cross section, and cross section spline (second derivative useful for spline interpolation).

Column	Type	Description
id	integer	Index (primary key)
element	text	Atomic symbol for element
log_energy	json array	log of Energy values (eV)
log_photoabsorption	json array	log of cross section ( $\text{cm}^2/\text{gr}$ )
log_photoabsorption_spline	json array	log of cross section spline

### 1.6.7 Scattering Table

The *scattering* table holds data for the coherent and incoherent X-ray scattering cross sections, in  $\text{cm}^2/\text{gr}$ . Each row represents an element.

DB Table of Coherent and Incoherent Scattering Cross Sections. JSON-encoded arrays are held for logs of energy, cross section, and cross section spline (second derivative useful for spline interpolation).

Column	Type	Description
id	integer	Index (primary key)
element	text	Atomic symbol for element
log_energy	json array	log of Energy values (eV)
log_coherent_scatter	json array	log of cross section ( $\text{cm}^2/\text{gr}$ )
log_coherent_scatter_spline	json array	log of cross section spline
log_incoherent_scatter	json array	log of cross section ( $\text{cm}^2/\text{gr}$ )
log_incoherent_scatter_spline	json array	log of cross section spline

### 1.6.8 Coster\_Kronig Table

The *Coster\_Kronig* table holds data for energy levels, partial and total transition probabilities for the Coster-Kronig transitions (Auger processes in which the empty core level is filled from an electron in a higher level with the same principle quantum number). The partial probability describes direct transitions, while the total probability includes cascade effects. Each row represents a transition.

DB Table of Coster-Kronig Transitions.

Column	Type	Description
id	integer	Index (primary key)
element	text	Atomic symbol for element
initial_level	text	IUPAC symbol for initial level
final_level	text	IUPAC symbol for final level
transition_probability	float	direct transition probability
total_transition_probability	float	total transition probability

### 1.6.9 Waasmaier Table

The *Waasmaier* table holds data for calculating elastic X-ray scattering factors  $f_0(k)$ , from [[Waasmaier and Kirfel (1995)]]. The scattering factor is unitless, and  $k = \sin(\theta)/\lambda$  where  $\theta$  is the scattering angle and  $\lambda$  is the X-ray wavelength. available for many common ionic states for each element. Each row represents an ion.

DB Table of Elastic Scattering Cross Section Coefficients

Column	Type	Description
id	integer	Index (primary key)
atomic_number	integer	Atomic Number, Z
element	text	Atomic symbol for element
ion	text	symbol for element and ionization
offset	float	offset value
scale	json array	coefficients for calculation
exponents	json array	coefficients for calculation

### 1.6.10 KeskiRahkonen\_Krause Table

The *KeskiRahkonen\_Krause* table holds data for energy widths of the core electronic levels from [[Keski-Rahkonen and Krause (1974)]]. Values are in eV, and each row represents an energy level for an element.

DB Table of Core Hole Widths from Keski-Rahkonen and Krause

Column	Type	Description
id	integer	Index (primary key)
atomic_number	integer	Atomic Number, Z
element	text	Atomic symbol for element
edge	text	IUPAC symbol for energy level ('K')
width	float	width of level (eV)

### 1.6.11 Krause\_Oliver Table

The *Krause\_Oliver* table holds data for energy widths of the core electronic levels from [\[Krause and Oliver \(1979\)\]](#). Values are in eV, and each row represents an energy level for an element.

DB Table of Core Hole Widths from Krause and Oliver

Column	Type	Description
id	integer	Index (primary key)
atomic_number	integer	Atomic Number, Z
element	text	Atomic symbol for element
edge	text	IUPAC symbol for energy level ('K')
width	float	width of level (eV)

### 1.6.12 Chantler Table

The *Chantler* table holds data for resonant X-ray scattering factors  $f'(E)$  and  $f''(E)$  as well as photo-electric absorption, coherent, and incoherent scattering factors from [\[Chantler \(2000\)\]](#). As with other tables, scattering factors are unitless, and cross sections are in cm<sup>2</sup>/gr. Each row represents an element.

DB Table of resonant scattering and mass attenuation coefficients from Chantler.

Column	Type	Description
id	integer	Index (primary key)
element	text	Atomic symbol for element
mue_f2	float	factor to convert $\mu(E)$ to $f''(E)$
density	float	atomic density (gr/cm <sup>3</sup> )
corr_henke	float	Henke correction to $f'(E)$
corr_cl35	float	Cromer-Liberman correction to $f'(E)$
corr_nucl	float	nuclear correction to $f'(E)$
energy	json array	energies for interpolation
f1	json array	$f'(E)$ (e)
f2	json array	$f''(E)$ (e)
mu_photo	json array	photoabsorption $\mu(E)$ (cm <sup>2</sup> /gr)
mu_incoh	json array	incoherent scattering (cm <sup>2</sup> /gr)
mu_total	json array	total attenuation (cm <sup>2</sup> /gr)

## 1.7 References



## BIBLIOGRAPHY

- [Als-Nielsen and McMorrow (2011)] J. Als-Nielsen and D. McMorrow. *Elements of Modern X-ray Physics, 2nd Edition*. John Wiley & Sons, 2011. URL: <https://dx.doi.org/10.1002/9781119998365>.
- [Chantler (2000)] C. T. Chantler. Detailed tabulation of atomic form factors, photoelectric absorption and scattering cross section, and mass attenuation coefficients in the vicinity of absorption edges in the soft X-ray ( $Z=30-36$ ,  $Z=60-89$ ,  $E=0.1$  keV-10 keV), addressing convergence issues of earlier work. *Journal of Physical and Chemical Reference Data*, 29(4):597–1048, JUL-AUG 2000. URL: <https://dx.doi.org/10.1063/1.1321055>.
- [Chantler (2016)] C. T. Chantler. FFAST Data on a finer energy grid: Personal Communication. 2016.
- [Elam, Ravel, and Sieber (2002)] W. T. Elam, B. D. Ravel, and J. R. Sieber. A new atomic database for x-ray spectroscopic calculations. *Radiation Physics and Chemistry*, 63(2):121–128, February 2002. URL: [https://dx.doi.org/10.1016/S0969-806X\(01\)00227-4](https://dx.doi.org/10.1016/S0969-806X(01)00227-4).
- [Hipp (2012)] D. R. Hipp. Sqlite. 2012. URL: <https://www.sqlite.org>.
- [Keski-Rahkonen and Krause (1974)] O. Keski-Rahkonen and M. O. Krause. Total and partial atomic-level widths. *Atomic Data and Nuclear Data Tables*, 14(2):139–146, 1974. URL: [https://dx.doi.org/10.1016/S0092-640X\(74\)80020-3](https://dx.doi.org/10.1016/S0092-640X(74)80020-3).
- [Knoll (2010)] G. F. Knoll. *Radiation Detection and Measurement, 4th Edition*. John Wiley & Sons, 2010. URL: <https://dx.doi.org/10.1002/9780470131480>.
- [Krause and Oliver (1979)] M. O. Krause and J. H. Oliver. Natural widths of atomic k and l levels, ka x-ray lines and several kll auger lines. *Journal of Physical and Chemical Reference Data*, 8:329, 1979. URL: <https://dx.doi.org/10.1063/1.555595>.
- [Waasmaier and Kirfel (1995)] D. Waasmaier and A. Kirfel. New analytical scattering factor functions for free atoms and ions. *Acta Crystallographica A*, 51:416–431, 1995. URL: <https://dx.doi.org/10.1107/S0108767394013292>.



## PYTHON MODULE INDEX

### X

xraydb, [16](#)





## A

`add_material()` (in module *xraydb*), 25  
`atomic_density()` (in module *xraydb*), 18  
`atomic_mass()` (in module *xraydb*), 18  
`atomic_number()` (in module *xraydb*), 18  
`atomic_symbol()` (in module *xraydb*), 18

## C

`chantler_energies()` (in module *xraydb*), 22  
`chemparse()` (in module *xraydb*), 24  
`ck_probability()` (in module *xraydb*), 21  
`coherent_cross_section_elam()` (in module *xraydb*), 22  
`core_width()` (in module *xraydb*), 19

## D

`darwin_width()` (in module *xraydb*), 27  
 DB Table of Basic Properties of the Elements, 36  
 DB Table of Coherent and Incoherent Scattering Cross Sections, 37  
 DB Table of Core Hole Widths, 38, 39  
 DB Table of Coster-Kronig Transitions, 38  
 DB Table of Database Versions, 36  
 DB Table of Elastic Scattering Cross Section Coefficients, 38  
 DB Table of Photoabsorption Cross Sections, 37  
 DB Table of resonant scattering and mass attenuation coefficients from Chantler, 39  
 DB Table of X-ray Levels, 36  
 DB Table of X-ray Transitions, 36

## F

`f0()` (in module *xraydb*), 18  
`f0_ions()` (in module *xraydb*), 18  
`f1_chantler()` (in module *xraydb*), 23  
`f2_chantler()` (in module *xraydb*), 23  
`find_material()` (in module *xraydb*), 25  
`fluor_yield()` (in module *xraydb*), 21

## G

`get_material()` (in module *xraydb*), 25  
`get_materials()` (in module *xraydb*), 24  
`get_xraydb()` (in module *xraydb*), 17  
`guess_edge()` (in module *xraydb*), 20

## I

`incoherent_cross_section_elam()` (in module *xraydb*), 22  
`ionchamber_fluxes()` (in module *xraydb*), 30  
`ionization_potential()` (in module *xraydb*), 29

## M

`material_mu()` (in module *xraydb*), 26  
`material_mu_components()` (in module *xraydb*), 27  
`mirror_reflectivity()` (in module *xraydb*), 29  
 module  
   *xraydb*, 16  
`mu_chantler()` (in module *xraydb*), 23  
`mu_elam()` (in module *xraydb*), 21

## T

Table of Effective Ionization Potentials, 10  
 Table of X-ray Edges, 33  
 Table of X-ray emission lines, 33

## V

`validate_formula()` (in module *xraydb*), 24

## X

X-ray Periodic Tables, 4  
`xray_delta_beta()` (in module *xraydb*), 27  
`xray_edge()` (in module *xraydb*), 19  
`xray_edges()` (in module *xraydb*), 19  
`xray_lines()` (in module *xraydb*), 20  
*xraydb*  
   module, 16  
*xraydb* Python module, 17